

# **MAD** *Analysis* **5**

## The LaTeX report

---

Generated by ernesto on 12 April 2019, 17:24:39

This report has been generated automatically by MADANALYSIS 5.

Please cite:

**E. Conte, B. Fuks and G. Serret,**  
*MadAnalysis 5, A User-Friendly Framework for Collider Phenomenology,*  
Comput. Phys. Commun. **184** (2013) 222-256,  
arXiv:1206.1599 [hep-ph].

To contact us:

<http://madanalysis.irmp.ucl.ac.be>  
ma5team@iphc.cnrs.fr

---

## Contents

<b>1</b>	<b>Setup</b>	<b>2</b>
1.1	Command history	2
1.2	Configuration	2
<b>2</b>	<b>Datasets</b>	<b>3</b>
2.1	dy	3
2.2	ma1000	3
2.3	dy_pt250	3
2.4	ma1000_pt250	4
<b>3</b>	<b>Histos and cuts</b>	<b>5</b>
3.1	Histogram 1	5
3.2	Histogram 2	6
3.3	Histogram 3	7

---

## 1 Setup

### 1.1 Command history

```
ma5>import /home/ernesto/Work/codes/MCgenerators/MadGraph/MG5_aMC_v2_updated/LFVHD/-  
DY/L0/Events/DY.lhco as DY  
ma5>import /home/ernesto/Work/codes/MCgenerators/MadGraph/MG5_aMC_v2_updated/LFVHD/-  
MSSM-Higgs2ditaL0/Events/MA1000.lhco as MA1000  
ma5>import /home/ernesto/Work/codes/MCgenerators/MadGraph/MG5_aMC_v2_updated/LFVHD/-  
DY/L0/Events/DY_pT250.lhco as DY_pT250  
ma5>import /home/ernesto/Work/codes/MCgenerators/MadGraph/MG5_aMC_v2_updated/LFVHD/-  
MSSM-Higgs2ditaL0/Events/MA1000_pT250.lhco as MA1000_pT250  
ma5>define tau = ta+ ta-  
ma5>define muon = mu+ mu-  
ma5>set main.stacking_method = superimpose  
ma5>set DY.linecolor = black  
ma5>set MA1000.linecolor = red  
ma5>set DY_pT250.linecolor = green  
ma5>set MA1000_pT250.linecolor = blue  
ma5>plot PT(tau) 50 0 500 [logY]  
ma5>plot M(tau) 50 0 150 [logY]  
ma5>plot M(tau) 50 0 20 [logY]  
ma5>submit Mtau_pTlep-cuts
```

### 1.2 Configuration

- MadAnalysis version 1.6.33 (2017/11/20).
- Histograms are not scaled.

## 2 Datasets

### 2.1 dy

- Samples stored in the directory: [/home/ernesto/Work/codes/MCgenerators/MadGraph/- MG5\\_aMC\\_v2\\_updated/HEPTools/madanalysis5/madanalysis5](#) .
- Sample consisting of: `signal` events.
- Generated events: `10000` events.
- Normalization to the luminosity: `0+- 0` events.
- Ratio (event weight): `0.0` .

Path to the event file	Nr. of events	Cross section (pb)	Negative wgts (%)
/home/ernesto/Work/codes/- MCgenerators/MadGraph/- MG5_aMC_v2_updated/- LFVHD/DY/LO/Events/DY.lhco	10000	0.0 @ 0.0%	0.0

### 2.2 ma1000

- Samples stored in the directory: [/home/ernesto/Work/codes/MCgenerators/MadGraph/- MG5\\_aMC\\_v2\\_updated/HEPTools/madanalysis5/madanalysis5](#) .
- Sample consisting of: `signal` events.
- Generated events: `10000` events.
- Normalization to the luminosity: `0+- 0` events.
- Ratio (event weight): `0.0` .

Path to the event file	Nr. of events	Cross section (pb)	Negative wgts (%)
/home/ernesto/Work/codes/- MCgenerators/MadGraph/- MG5_aMC_v2_updated/- LFVHD/MSSM-Higgs2ditaui/- LO/Events/MA1000.lhco	10000	0.0 @ 0.0%	0.0

### 2.3 dy\_pt250

- Samples stored in the directory: [/home/ernesto/Work/codes/MCgenerators/MadGraph/- MG5\\_aMC\\_v2\\_updated/HEPTools/madanalysis5/madanalysis5](#) .
- Sample consisting of: `signal` events.
- Generated events: `10000` events.
- Normalization to the luminosity: `0+- 0` events.

- Ratio (event weight): 0.0 .

Path to the event file	Nr. of events	Cross section (pb)	Negative wgts (%)
/home/ernesto/Work/codes/-MCgenerators/MadGraph/-MG5_aMC_v2_updated/-LFVHD/DY/LO/Events/-DY_pT250.lhco	10000	0.0 @ 0.0%	0.0

## 2.4 ma1000\_pt250

- Samples stored in the directory: /home/ernesto/Work/codes/MCgenerators/MadGraph/-MG5\_aMC\_v2\_updated/HEPTools/madanalysis5/madanalysis5 .
- Sample consisting of: signal events.
- Generated events: 10000 events.
- Normalization to the luminosity: 0+- 0 events.
- Ratio (event weight): 0.0 .

Path to the event file	Nr. of events	Cross section (pb)	Negative wgts (%)
/home/ernesto/Work/codes/-MCgenerators/MadGraph/-MG5_aMC_v2_updated/-LFVHD/MSSM-Higgs2ditaui/-LO/Events/MA1000_pT250.lhco	10000	0.0 @ 0.0%	0.0

### 3 Histos and cuts

#### 3.1 Histogram 1

\* Plot: PT ( tau )

Dataset	Integral	Entries per event	Mean	RMS	% underflow	% overflow
dy	4191	1.15	37.8416	12.02	0.0	0.0
ma1000	8043	1.26	269.837	125.5	0.0	2.337
dy_pt250	7788	1.25	220.677	107.1	0.0	1.746
ma1000_pt250	7986	1.27	281.549	120.8	0.0	2.504




Figure 1.

### 3.2 Histogram 2

\* Plot: M ( tau )

Dataset	Integral	Entries per event	Mean	RMS	% underflow	% overflow
dy	4191	1.15	3.47246	2.648	1.623	0.0
ma1000	8043	1.26	11.9557	12.79	3.419	0.0
dy_pt250	7788	1.25	10.4759	10.34	3.017	0.0
ma1000_pt250	7986	1.27	12.3108	12.69	3.306	0.0




Figure 2.

### 3.3 Histogram 3

\* Plot: M ( tau )

Dataset	Integral	Entries per event	Mean	RMS	% underflow	% overflow
dy	4191	1.15	3.47246	2.648	1.623	0.02386
ma1000	8043	1.26	11.9557	12.79	3.419	16.28
dy_pt250	7788	1.25	10.4759	10.34	3.017	11.68
ma1000_pt250	7986	1.27	12.3108	12.69	3.306	17.56




Figure 3.